



Recombinant *Saccharomyces cerevisiae*

Uncharacterized protein YJL213W (YJL213W)

Product Code	CSB-MP336800SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40896
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MLMPKSVALL RIRYTLAAL ARGFTTVRDC GGAEGFLKAE IRQGSLLNGPR LITCGHAISQ TGGHGDLSRG ALPASAFDSC SCHFGQVGVV ADGVPECYKA AREEFRRGAD FIKIMGGGGV ASPTDKISNK QFCDDEIKAL VDVANSYHTY VTAHAYTAEA IQNCIKLGVK GIEHGNLLDE RTAELMAELG CYLTPTLVTY KVMGSDQFSA FLGPENSRKN TEVLYKGIDA MKIAQKKKVK ICFGSDLLGP LYGYQTQEFR IRGKVQTTQE ILLSATVTPA EMNGLGDKLG QIKPGFIADL LMMKSNPLDD ITILDEPESN ILFVMKEGRI Y
Source	Mammalian cell
Target Names	YJL213W
Protein Names	Recommended name: Uncharacterized protein YJL213W
Expression Region	1-331
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.