



Recombinant *Saccharomyces cerevisiae* AP-1-like transcription factor YAP4 (CIN5)

Product Code	CSB-YP336802SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P40917
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MLMQIKMDNH PFNFQPILAS HSMTRDSTKP KKMTDTAFVP SPPVGFIIKEE NKADLHTISV VASNVTLPQI QLPKIATLEE PGYESRTGSL TDLSGRRNSV NIGALCEDVP NTAGPHIARP VTINNLIPPS LPRLNTYQLR PQLSDTHLNC HFNSNPYTTA SHAPFESSYT TASTFTSQPA ASYFPSNSTP ATRKNSATTN LPSEERRRVS VSLSEQVFNE GERYNNDGQL IGKTGKPLRN TKRAAQNRSA QKAQRQREK YIKNLEEKSK LFDGLMKENS ELKKMIESLK SKLKE
Source	Yeast
Target Names	CIN5
Protein Names	Recommended name: AP-1-like transcription factor YAP4 Alternative name(s): Chromosome instability protein 5 Transcription activator CIN5
Expression Region	1-295
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.