



Recombinant *Saccharomyces cerevisiae* Glycerol-3-phosphate dehydrogenase [NAD (+)] 2, mitochondrial

Product Code	CSB-BP336933SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P41911
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	HPVL YTRRAYKILP SRSTFLRRSL LQTQLHSKMT AHTNIKQHKH CHEDHPIRRS DSAVSIVHLK RAPFKVTVIG SGNWGTIAK VIAENTELHS HIFEPEVRMW VFDEKIGDEN LTDIINTRHQ NVKYLPNIDL PHNLVADPDL LHSIKGADIL VFNIPHQFLP NIVKQLQGHV APHVRAISCL KGFELGSKGV QLLSSYVTDE LGIQCGALSG ANLAPEVAKE HWSETTVAYQ LPKDYQGDGK DVDHKILKLL FHRPYFHVNV IDDVAGISIA GALKNVVALA CGFVEGMGWG NNASAAIQRL GLGEIKFGR MFFPESKVET YYQESAGVAD LITTCGGRN VKVATYMAKT GKSALEAEKE LLNGQSAQGI ITCREVHEWL QTCELTQEFP LFEAVYQIVY NNVRMEDLPE MIEELDIDDE
Source	Baculovirus
Target Names	GPD2
Protein Names	Recommended name: Glycerol-3-phosphate dehydrogenase [NAD(+)] 2, mitochondrial EC= 1.1.1.8
Expression Region	17-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.