



Recombinant Solanum tuberosum Mitochondrial outer membrane protein porin of 36 kDa

Product Code	CSB-BP336972FIG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P42056
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	VKGPGLYSD IGKKARDLLY RDYVSDHKFT VTTYSTTGVA ITASGLKKGE LFLADVSTQL KNKNITTDVK VDTNSNVYTT ITVDEPAPGL KTIFSFVVPD QKSGKVELQY LHEYAGINTS IGLTASPLVN FSGVAGNNTV ALGTDLSFDT ATGNFTKCNA GLSFSSSDLI ASLALNDKGD TVSASYHTV KPVTNTAVGA ELTHSFSSNE NTLTIGTQHL LDPLTTVKAR VNSYGKASAL IQHEWRPKSL FTISGEVDTR AIEKSAKIGL AVALKP
Source	Baculovirus
Protein Names	Recommended name: Mitochondrial outer membrane protein porin of 36 kDa Alternative name(s): POM 36 Voltage-dependent anion-selective channel protein Short name= VDAC
Expression Region	2-276
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.