



# Recombinant Solanum tuberosum Mitochondrial outer membrane protein porin of 36 kDa

<b>Product Code</b>	CSB-YP336972FIG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P42056
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Solanum tuberosum (Potato)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	VKGPGLYSD IGKKARDLLY RDYVSDHKFT VTTYSTTGVA ITASGLKKGE LFLADVSTQL KNKNITTDVK VDTNSNVYTT ITVDEPAPGL KTIFSFVVPD QKSGKVELQY LHEYAGINTS IGLTASPLVN FSGVAGNNTV ALGTDLSFDT ATGNFTKCNA GLSFSSDLI ASLALNDKGD TVSASYHTV KPVTNTAVGA ELTHSFSSNE NTLTIGTQHL LDPLTTVKAR VNSYGKASAL IQHEWRPKSL FTISGEVDTR AIEKSAKIGL AVALKP
<b>Source</b>	Yeast
<b>Protein Names</b>	Recommended name: Mitochondrial outer membrane protein porin of 36 kDa Alternative name(s): POM 36 Voltage-dependent anion-selective channel protein Short name= VDAC
<b>Expression Region</b>	2-276
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.