



Recombinant Ribonucleoside-diphosphate reductase small chain (rnr-2)

Product Code	CSB-BP336991CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P42170
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MTLTEIQNVE KENAGASVPK HSSNKLKLEK ELEKLEIVDQ TKAASAEETN NESEVNELDA DEPMLQDLDN RRVIFPLKHH DIWNFYKKAV ASFWTVVEVD LGKDMNDWEK MNGDEQYFIS RILAFFAASD GIVNENLCER FSNEVQVSEA RFFYGFQIAI ENIHSEMYSK LIETYIRDET ERNTLFNAVD EFEFIKKKAD WALRWISDKK ASFAERLIAF AAVEGIFFSG SFASIFWLKK RGLMPGLTHS NELISRDEGL HRDFACLLYS KLQKKLQQR IYDIKDAVA IEQEFLTEAL PVDMIGMNCR LMSQYIEFVA DHLLVELGCD KLYKSKNPFDFMENISIDGK TNFFEKRVSE YQRPQVMVNE AERQFDLEAD F
Source	Baculovirus
Target Names	rnr-2
Protein Names	Recommended name: Ribonucleoside-diphosphate reductase small chain EC=1.17.4.1 Alternative name(s): Ribonucleotide reductase small subunit
Expression Region	1-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.