



Recombinant Manganese ABC transporter substrate-binding lipoprotein (psaA)

Product Code	CSB-YP337054FMW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P42363
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pneumoniae
Purity	>85% (SDS-PAGE)
Sequence	CSKGTSSKSS SDKLKVVTTN SILADITKNI AGDKIELHSI VPVGQDPHEY EPLPEDVKKT SQADLIFYNG INLETGGNAW FTKLVKNANK VENKDYFAAS DGVEVIYLEG QNQAGKEDPH AWLNLENGII YAKNIAKQLI AKDPKNKDFY EKNLAAYTEK LSKLDQEAKQ AFNNIPAEEK MIVTSEGCFK YFSKAYGVPS AYIWEINTEV EGTPEQIKTL LEKLRQTKVP SLFVESSVDE RPMKTVSKDS NIPIFAKIFT DSIAKEGEEG DSYYSMMKWN LEKIAEGLNK
Source	Yeast
Target Names	psaA
Protein Names	Recommended name: Manganese ABC transporter substrate-binding lipoprotein Alternative name(s): Pneumococcal surface adhesin A
Expression Region	21-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.