



Recombinant Arabidopsis thaliana Ethylene-responsive transcription factor RAP2-3 (RAP2-3)

Product Code	CSB-MP337192DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P42736
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MCGGAIISDY APLVTKAKGR KLTAEELWSE LDASAADDFW GFYSTSKLHP TNQVNVKEEA VKKEQATEPG KRRKRKNVYR GIRKRPWGKW AAEIRDPRKG VRVWLGTFTNT AEEAAMAYDV AAKQIRGDKA KLNFPDLHHP PPPNYTPPPS SPRSTDQPPA KKVCVVSQSE SELSQPSFPV ECIGFGNGDE FQNLSYGFEP DYDLKQQISS LESFLELDGN TAEQPSQLDE SVSEVDMWML DDVIASYE
Source	Mammalian cell
Target Names	RAP2-3
Protein Names	Recommended name: Ethylene-responsive transcription factor RAP2-3 Alternative name(s): Cadmium-induced protein AS30 Ethylene response factor 72 Short name= ERF72 Ethylene-responsive element binding protein Short name= AtEBP
Expression Region	1-248
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.