



Recombinant Arabidopsis thaliana Glutathione S-transferase F6 (GSTF6)

| | |
|--------------------------|---|
| Product Code | CSB-YP337200DOA |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P42760 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MAGIKVFGHP ASTATRRVLI ALHEKNVDFE FVHVELKDGE HKKEPFILRN PFGKVPAFED GDFKIFESRA ITQYIAHEFS DKGNNLLSTG KDMAIIAMGI EIESHEFDPV GSKLVWEQVL KPLYGMTTDDK TVVEEEEAKL AKVLDVYEHR LGESKYLASD HFTLVDLHTI PVIQYLLGTP TKKLFDERPH VSAWVADITS RPSAQKVL |
| Source | Yeast |
| Target Names | GSTF6 |
| Protein Names | Recommended name: Glutathione S-transferase F6 Short name= AtGSTF6 EC= 2.5.1.18 Alternative name(s): AtGSTF3 GST class-phi member 6 Glutathione S-transferase 1 Short name= AtGST1 Protein EARLY RESPONSE TO DEHYDRATI |
| Expression Region | 1-208 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |