



Recombinant Escherichia coli Putative tagatose-6-phosphate ketose/aldose isomerase (agaS)

Product Code	CSB-YP337256ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P42907
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MPENYTPAAA ATGTWTEEEI RHQPRAWIRS LTNIDALRSA LNNFLEPLLK KENLRIILTG AGTSAFIGDI IAPWLASHTG KNFSAVPTTD LVTNPMDYLN PAHPLLLISF GRSGNSPESV AAVELANQFV PECYHLPITC NEAGALYQNA INSDNAFALL MPAETHDRGF AMTSSITMM ASCLAVFAPE TINSQTFRDV ADRCQAILTS LGDFSEGVFG YAPWKRIVYL GSGGLQGAAR ESALKVLELT AGKLAIFYDS PTGFRHGPKS LVDDETLVVV FVSSHPYTRQ YDLDLLAELR RDNQAMRVIA IAAESSDIVA AGPHIILPPS RHFIDVEQAF CFLMYAQTF LMQSLHMGNT PDTPSASGTV NRVVQGVIIH PWQA
Source	Yeast
Target Names	agaS
Protein Names	Recommended name: Putative tagatose-6-phosphate ketose/aldose isomerase EC= 5.3.1.-
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.