



Recombinant Phleum pratense Pollen allergen Phl p 2 (PHLP2)

Product Code	CSB-YP337346EUQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P43214
Product Type	Recombinant Protein
Immunogen Species	Phleum pratense (Common timothy)
Purity	>85% (SDS-PAGE)
Sequence	VPKV TFTVEKGSNE KHLAVLVKYE GDTMAEVELR EHGSDWVAM TKGEGGVWTF DSEEPLQGP F NFRFLTEKGM KNVFDDVPE KYTIGATYAP EE
Source	Yeast
Target Names	PHLP2
Protein Names	Recommended name: Pollen allergen Phl p 2 Alternative name(s): Allergen Phl p II Allergen= Phl p 2
Expression Region	27-122
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.