



Recombinant Arabidopsis thaliana Transcription factor TGA2 (TGA2)

Product Code	CSB-EP337365DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P43273
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MADTSPRTDV STDDTDHPD LGSEGALVNT AASDSSDRSK GKMDQKTLRR LAQNREAARK SRLRKKAYVQ QLENSRLKLT QLEQELQRAR QQGVFISGTG DQAHSTGGNG ALAFDAEHSR WLEEKKNQMN ELRSALNAHA GDSELRIVD GVMAHYEELF RIKSNAAKND VFHLLSGMWK TPAERCFLWL GGFRSSELLK LLANQLEPMT ERQLMGINNLL QQTSQQAEDA LSQGMESLQQ SLADTLSSGT LGSSSSGNVA SYMGQMAMAM GKLGTLEGFI RQADNLRQLT LQQMIRVLTT RQSARALLAI HDYFSRLRAL SSLWLARPRE
Source	E.coli
Target Names	TGA2
Protein Names	Recommended name: Transcription factor TGA2 Alternative name(s): HBP-1b homolog Short name= AHBP-1b bZIP transcription factor 20 Short name= AtbZIP20
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.