



Recombinant *Saccharomyces cerevisiae* LAS seventeen-binding protein 3 (LSB3)

Product Code	CSB-EP337461SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P43603
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	<p> MGINNPIPRS LKSETKKA AK ILASFVKPNQ VFGADQVIPP DVLKRAKGLA IITILKAGFL FSGRAGSGVI VARLKDGTWS APSAIAMAGA GAGGMVGIEL TDFVFILNTQ DAVKSFSEFG TITLGGNVSV SAGPLGRSAE AAASASAGGV AAVFAYS KSK GLFAGV SVEG SAIERREAN RKFYGDNCTA KMILSGRIRP PPAVDPLFRV LESRAFNYRP SNGGRGSFDD DEDDYDDDD YYNDIPSSFS STDASSTRPN TRSTRRAQS GSRYTFDDDD DDDDYGTGYS RNSRLAPTNS GGSGGKLDDP SGASSYYASH RRS GTAQSRA RSSRNWADD EYDDYDDDYE SGYRRGN GRD RTKDREVDL SNRFSKSRIS SASTPQTSQG RFTAPTSPST SSPKAVALYS FAGEESGDL P FRKGDVITIL KKSDSQNDWW TGRVNGREGI FPANYVELV </p>
Source	E.coli
Target Names	LSB3
Protein Names	Recommended name: LAS seventeen-binding protein 3 Short name= LAS17-binding protein 3
Expression Region	1-459
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.