



Recombinant *Thermus thermophilus* Uridylate kinase (pyrH)

Product Code	CSB-EP337567TNT-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P43891
Product Type	Recombinant Protein
Immunogen Species	<i>Thermus thermophilus</i> (strain HB8 / ATCC 27634 / DSM 579)
Purity	>85% (SDS-PAGE)
Sequence	MKYKRVLLKL SGEFLTRNGF GIEPEATQAL AREIKAAAYDT GVQLAIVIGA GNLWRGARQG VGMDRATADY IGMLATIMNA LALQDALESL GVPTRVQTAL TITQVAEPYI RRRALRHLEK ERIVIFGGGT GNPFFSTDTA AALRALEVGA EVLMAKNKV DGVYSDDPRK NPEAVRFDEL TYLEVLNRGL QVMDTTAITL CMEAGLPIVV FDIFKPGALV GIIQGEKVGTL LIH
Source	E.coli
Target Names	pyrH
Protein Names	Recommended name: Uridylate kinase Short name= UK EC= 2.7.4.22 Alternative name(s): Uridine monophosphate kinase Short name= UMP kinase Short name= UMPK
Expression Region	1-233
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.