




# Recombinant Human Butyrophilin subfamily 3 member A1 (BTN3A1), partial

<b>Product Code</b>	CSB-YP002873HU
<b>Relevance</b>	Plays a role in T-cell activation and in the adaptive immune response. Regulates the proliferation of activated T-cells. Regulates the release of cytokines and IFNG by activated T-cells. Mediates the response of T-cells toward infected and transformed cells that are characterized by high levels of phosphorylated metabolites, such as isopentenyl pyrophosphate.
<b>Abbreviation</b>	Recombinant Human BTN3A1 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O00481
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	QFSVLGSPGPILAMVGEDADLPCHLFPTMSAETMELKWWSSSLRQVVNVYAD GKEVEDRQSAPYRGRTSILRDGITAGKAALRIHNVTASDSGKYLQYFQDGDY EKALVELKVAALGSDLHVDVKGYKDGGIHLECRSTGWYPQPQIQWSNNKGEN IPTVEAPVVADGVGLYAVAASVIMRGSSGEGVSVCTIRSSLLGLEKTASISADPF FRSAQRWIAALAG
<b>Research Area</b>	Immunology
<b>Source</b>	Yeast
<b>Target Names</b>	BTN3A1
<b>Protein Names</b>	CD277
<b>Expression Region</b>	30-254aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	26.2 kDa
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.