



Recombinant Haemophilus influenzae FKBP-type peptidyl-prolyl cis-trans isomerase slyD (slyD)

Product Code	CSB-YP337956HTA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P44830
Product Type	Recombinant Protein
Immunogen Species	Haemophilus influenzae (strain ATCC 51907 / DSM 11121 / KW20 / Rd)
Purity	>85% (SDS-PAGE)
Sequence	MKVEKNVVVS ISYQVRTQDG VLVDEAPANQ PLEYLQGHNN LVIGLEKALE GKEVGDKFEV RVQPEEGYGA YSENMVQRVP KDVFGQVDEL EVGMRFLADT DIGPVPVVIT EIDGDEVVVD GNHMLAGQEL HFTVEVVAAR EATLEEIAHG HVHGAHSHDD DEEGHGCGCG GHHHEHNHEH NHGSCGCGGH
Source	Yeast
Target Names	slyD
Protein Names	Recommended name: FKBP-type peptidyl-prolyl cis-trans isomerase slyD Short name= PPIase EC= 5.2.1.8 Alternative name(s): Metallochaperone slyD
Expression Region	1-190
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.