



# Recombinant Human CCN family member 1 (CCN1)

<b>Product Code</b>	CSB-BP006463HU
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O00622
<b>Storage Buffer</b>	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	TCPAAC HCPLEAPKCA PGVGLVRDGC GCCKVCAKQL NEDCSKTQPC DHTKGLECNF GASSTALKGI CRAQSEGRPC EYNSRIYQNG ESFQPNCXHQ CTCIDGAVGC IPLCPQELSL PNLGCPNPR L VKVTGQCCEE WVCDEDSIKD PMEDQDGLLG KELGFDASEV ELTRNNELIA VGKGSLLKRL PVFGMEPRIL YNPLQGQKCI VQTTSWSQCS KTCGTGISTR VTNDNPECRL VKETRICEVR PCGQPVYSSL KKGKKCSKTK KSPEPVRFTY AGCLSVKKYR PKYCGSCVDG RCCTPQLTRT VKMRFRCEDG ETFSKNVMMI QSCKCNYNCP HANEAAPPFY RLFNDIHKFR D
<b>Source</b>	Baculovirus
<b>Target Names</b>	CYR61
<b>Protein Names</b>	Recommended name: Protein CYR61 Alternative name(s): CCN family member 1 Cysteine-rich angiogenic inducer 61 Insulin-like growth factor-binding protein 10 Short name= IBP-10 Short name= IGF-binding protein 10 Short name=
<b>Expression Region</b>	25-381
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.