



Recombinant Human CCN family member 1 (CCN1)

Product Code	CSB-MP006463HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O00622
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	TCPAAC HCPLEAPKCA PGVGLVRDGC GCCKVCAKQL NEDCSKTQPC DHTKGLECNF GASSTALKGI CRAQSEGRPC EYNSRIYQNG ESFQPNCXHQ CTCIDGAVGC IPLCPQELSL PNLGCPNPR L VKVTGQCCEE WVCDEDSIKD PMEDQDGLLG KELGFDASEV ELTRNNELIA VGKGSLLKRL PVFGMEPRIL YNPLQGQKCI VQTTSWSQCS KTCGTGISTR VTNDNPECRL VKETRICEVR PCGQPVYSSL KKGKKCSKTK KSPEPVRFTY AGCLSVKKYR PKYCGSCVDG RCCTPQLTRT VKMRFRCEDG ETFSKNVMMI QSCKCNYNCP HANEAAPPFY RLFNDIHKFR D
Source	Mammalian cell
Target Names	CYR61
Protein Names	Recommended name: Protein CYR61 Alternative name(s): CCN family member 1 Cysteine-rich angiogenic inducer 61 Insulin-like growth factor-binding protein 10 Short name= IBP-10 Short name= IGF-binding protein 10 Short name=
Expression Region	25-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.