



Recombinant Human Netrin-3 (NTN3)

Product Code	CSB-MP016128HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O00634
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	DPC HDEGGAPRGC VPGLVNAALG REVLASSTCG RPATRACDAS DPRRAHSPAL LTSPGGTASP LCWRSESLPR APLNVTLTVP LGKAFELVFV SLRFCSAPPA SVALLKSQDH GRSWAPLGFF SSHCDLDYGR LPAPANGPAG PGPEALCFPA PLAQPDGSL LAFSMQDSSP PGLDLDSSPV LQDWVTATDV RVVLTRPSTA GPRDMEAVV PYSYAATDLQ VGGRCCKNGH ASRCLLDTQG HLICDCRHGT EGPDCGRCKP FYCDRPWQRA TARESHACLA CSCNGHARRC RFNMELYRLS GRRSGGVCLN CRHNTAGRHC HYCREGFYRD PGRALSDRRA CRACDCHPVG AAGKTCNQTT GQCPCCKDGVV GLTCNRCAPG FQQRSPVAP CVKTIPIGPPT EDSSPVQPQD CDSHCKPARG SYRISLKKFC KKDYAVQVAV GARGEARGAW TRFPVAVLAV FRSGEERARR GSSALWVPAG DAACGCPRLP PGRRYLLGG GPGAAAGGAG GRGPGLIAAR GSLVLPWRDA WTRRLRRLQR RERRGRCSAA
Source	Mammalian cell
Target Names	NTN3
Protein Names	Recommended name: Netrin-3 Alternative name(s): Netrin-2-like protein
Expression Region	28-580
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.