



Recombinant Mouse Ephrin-A3 (Efna3)

Product Code	CSB-YP007462MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O08545
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	QGGPGALG NRHAVYWNSS NQHLLRREGYT VQVNVNDYLD IYCPHYNSSG PGGGAEQYVL YMVNLSGYRT CNASQGSKRW ECNRQHASHS PIKFSEKFQR YSAFSLGYEF HAGQEYYYIS TPTHNLHWKC LRMKVFVCCA STSHSGEKPV PTLPQFTMGP NVKINVLEDF EGENPQVPKL EKSISG
Source	Yeast
Target Names	Efna3
Protein Names	Recommended name: Ephrin-A3 Alternative name(s): EHK1 ligand Short name= EHK1-L EPH-related receptor tyrosine kinase ligand 3 Short name= LERK-3
Expression Region	23-206
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a member of the ephrin (EPH) family. The ephrins and EPH-related receptors comprise the largest subfamily of receptor protein-tyrosine kinases and have been implicated in mediating developmental events, especially in the nervous system and in erythropoiesis. Based on their structures and sequence relationships, ephrins are divided into the ephrin-A (EFNA) class, which are anchored to the membrane by a glycosylphosphatidylinositol linkage, and the ephrin-B (EFNB) class, which are transmembrane proteins. This gene encodes an EFNA class ephrin.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.