



Recombinant *Campylobacter jejuni* subsp. *jejuni* serotype O:2 Uncharacterized protein Cj0990c (Cj0990c)

Product Code	CSB-YP338202CZA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P45489
Product Type	Recombinant Protein
Immunogen Species	<i>Campylobacter jejuni</i> subsp. <i>jejuni</i> serotype O:2 (strain ATCC 700819 / NCTC 11168)
Purity	>85% (SDS-PAGE)
Sequence	MMMKIVFKGK SSEYEIQRSC FCVDAFVIKD KIEERDGVDF ITSNVDLLEF SDDSFTFEEI VKHFNICDTE DMIIVEDFDM KSNKDNQSQE DDIEHNILKS EKIIHENTKQ TSMQFKNLKF FSRIFKNENF LSDFKESKQE VVTIKKHEKL EIFKNLSQED QEISFVKIEI LNYDSNEDSL SFNLDIFPSG MSYKYGILKG SMHILQGKT SSTMLFPFLK SMIYKNKSEN SSEKIFTLMI NQKKHYKLIA NLS
Source	Yeast
Target Names	Cj0990c
Protein Names	Recommended name: Uncharacterized protein Cj0990c
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.