



Recombinant *Campylobacter jejuni* subsp. *jejuni* serotype O:2 Uncharacterized lipoprotein Cj0983 (Cj0983)

Product Code	CSB-EP338203CZA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P45492
Product Type	Recombinant Protein
Immunogen Species	<i>Campylobacter jejuni</i> subsp. <i>jejuni</i> serotype O:2 (strain ATCC 700819 / NCTC 11168)
Purity	>85% (SDS-PAGE)
Sequence	CGN SIDEKTVKKY ENQLNQTVKQ EIASLSQDSG IKIEFSDFKC NADGDFIACL SPNFKLTKD NNDEYQELFQ AKNIKIRSNE IYKGGANASI SIKEYYNDLF KNQKSIQSNL VFEDFKLGEK VVSDINASLF QQDPKISSFI NKLSSDSYTL SFDNSINKQE NNYLDNLDIK FYNAKLNFNNT NLNINLKEDL LNYLDSKGIK FNTQTLVMDK QAINELLNIA NYEQASDFSN TIQKYIILNN FKIDSTLKTE GVFSSYIATA KENLQTLKAQ SQNEEQALIF DKALAILNNI TQNDDYKLNLDLKFKNIPVS DYSTQGIDSI EKLSINNQDA TEALKIILPF IMFSMLMGGA SF
Source	E.coli
Target Names	Cj0983
Protein Names	Recommended name: Uncharacterized lipoprotein Cj0983
Expression Region	18-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.