



Recombinant Escherichia coli Uncharacterized protein yfeN (yfeN)

Product Code	CSB-BP338226ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P45564
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MKKHLLTLTL SSILAI PVVS HAEFKGGFAD IGVHYLDWTS RTTEKSSTKS HKDDFGYLEF EGGANFSWGE MYGFFDWENF YNGRHNKPGS EQRYTFKNTN RIYLGDTGFN LYLHAYGTYG SANRVNFHDD MFLYGIGYNF TGSGWWFKPF FAKRYTDQTY YTGDNQYVAG WVAGYNFMLG SEKFTLTNWN EYEFDRDATY AAGNGGKEGL NGAVALWWNA TSHITTGIQY RYADDKLGED FYQDAIYSI KFNF
Source	Baculovirus
Target Names	yfeN
Protein Names	Recommended name: Uncharacterized protein yfeN
Expression Region	1-254
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.