



# Recombinant Mouse Recombining binding protein suppressor of hairless-like protein (Rbpjl)

<b>Product Code</b>	CSB-EP019487MO-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O08674
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDPRETTDPS LPPGPLTHLS LPDSSEVRLQ SDGPSLLGSW TRSPPEHAI LREGVRTCLQ QRCEQTVWIL HAKVAQKSYG NEKRFFCPPP CVYLAGPGWR VKPMQDQALQ SAETGPTVCG YMGLDGASGS APETQKLNFE EQPDSREFGC AKTLYISDAD KRKHFRLLVLR LVLRRGGQELG TFHSRLIKVI SKPSQKKQSL KNTDLCISSG SKVSLFNRLR SQTVSTRYLS VEDGAFVASA RQWAAFTLHL ADGHCSQGDF PPQEGYIRYG SLVQLVCTVT GITLPPMIIR KVAKQCALLD VDEPISQLHK CAFQFPSDTP GGAGTYLCLA TEKVRFQAS LCPKEANRAL LNDSSCWTTI GTESVEFTFS TSLACTREPV TPVPLISTLE LSGGGDVATL ELHGENFHAG LKVWFGDVEA ETMYRSPRSL VCVVPDVAAF GSDWRWLRTP ITVPVLLRA DALFYSPFSS FTYTPEYSAL PRLPNAQEPA PDADTLLESI HHEFTRTNFH LFCPT
<b>Source</b>	E.coli
<b>Target Names</b>	Rbpjl
<b>Protein Names</b>	Recommended name: Recombining binding protein suppressor of hairless-like protein Alternative name(s): Transcription factor RBP-L
<b>Expression Region</b>	1-515
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	In mouse, recombining binding protein L (RBP-L) is a transcription factor that binds to DNA sequences almost identical to that bound by the Notch receptor signalling pathway transcription factor RBP-J. However, unlike RBP-J, RBP-L does not interact with Notch receptors. RBP-L has been shown to activate transcription in concert with Epstein-Barr virus nuclear antigen-2 (EBNA2). This protein is similar in sequence to the mouse RBP-L protein and Drosophila suppressor of hairless protein.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.