



Recombinant *Rhodothermus marinus* Beta-glucanase (bglA)

| | |
|--------------------------|---|
| Product Code | CSB-EP338318RMC |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P45798 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Rhodothermus marinus</i> (<i>Rhodothermus obamensis</i>) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | SDRSDKAPHW ELVWSDEFDY SGLPDPEKWD YDVGGHGWGN QELQYYTRAR IENARVGGGV LIIEARHEPY EGREYTSARL VTRGKASWTY GRFEIRARLP SGRGTWPAIW MLPDRQTYGS AYWPDNGEID IMEHVGFNPD VVHGTVHTKA YNHLLGTQRG GSIRVPTART DFHVYAIEWT PEEIRWFVDD SLYYRFPNER LTDPEADWRH WPFDQPFHLI MNIAVGGAWG GQQGVDPEAF PAQLVVDYVR VYRWVE |
| Source | E.coli |
| Target Names | bglA |
| Protein Names | Recommended name: Beta-glucanase EC= 3.2.1.73 Alternative name(s): 1,3-1,4-beta-D-glucan 4-glucanohydrolase Endo-beta-1,3-1,4 glucanase Lichenase |
| Expression Region | 31-286 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |