



Recombinant Escherichia coli HTH-type transcriptional regulator dsdC (dsdC)

Product Code	CSB-BP338417ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46068
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MEPLREIRNR LLNGWQLSKM HTFEVAARHQ SFALAAEELS LSPSAVSHRI NQLEELGIQ LFVRSHRKVE LTHEGKRVYW ALKSSLDTLN QEILDIKNQE LSGTLTLYSR PSIAQCWLVP ALGDFTRRYP SISLTVLTGN DNVNLQRAGI DLAIYFDDAP SAQLTHHFLM DEEILPVCSP EYAQRHALTN TVINLCHCTL LHDRQAWSND SGTDEWHSWA QHYAVNLPTS SGIGFDRSDL AVIAAMNHIG VAMGRKRLVQ KRLASGELVA PFGDMTVKCH QHYIYITLPG RQWPKEAFI IWLREQVKT S
Source	Baculovirus
Target Names	dsdC
Protein Names	Recommended name: HTH-type transcriptional regulator dsdC Alternative name(s): D-serine deaminase activator
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.