



Recombinant Escherichia coli Putative acyl-CoA thioester hydrolase ybhC (ybhC)

Product Code	CSB-YP338436ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46130
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	CSSTPPDQR PSDQTAPGTS SRPILSAKEA QNFDAQHYFA SLTPGAAAWN PSPITLPAQP DFVVGPAQTQ GVTHHTTIQAA VDAAIKRTN KRQYIAVMGP EYQGTVYVPA APGGITLYGT GEKPIDVKIG LSLDGGMSPA DWRHDVNPRG KYMPGKPAWY MYDSCQSKRS DSIGVLC SAV FWSQNNGLQL QNLTIENTLG DSVDAGNHPA VALRTDGDQV QINNVNIGR QNTFFVTNSG VQNRLETNRQ PRTLVTNSYI EGDVDIVSGR GAVVFDNTEF RVVNSRTQQE AYVFAPATLS NIYYGFLAVN SRFNAFGDV AQLGRSLDVD ANTNGQVVIR DSAINEGFNT AKPWADAVIS NRPFAGNTGS VDDNDEIQRN LNDTNYNRMW EYNNRGVGSK VVAEAKK
Source	Yeast
Target Names	ybhC
Protein Names	Recommended name: Putative acyl-CoA thioester hydrolase ybhC EC= 3.1.2.-
Expression Region	22-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.