



Recombinant Mouse SH2B adapter protein 3 (Sh2b3)

Product Code	CSB-MP021208MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O09039
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MNEPTVQPSR TSSAPASPAS PRGWSDFCEQ HAAAAARELA RQYWLFAH PQQPRADLVS LQFAELFQRH FCREVRESLA GPPGHDIYRAT APPRPALPKA RSSEDLGPRP ACALQHLRRG LRQLFRRRSA GELPGATSDT NDIDTTAASR PGPARKLLPW GLREPTEAL KEVVLRYSLA DEAAMDSGAR WQRGRLVLR PGPGHSHFLQ LFDPPKSSKP KLQEACSSIR EVRPCTRLEM PDNLYTFVLK VQDQTDIIFE VGDEQQLNSW LAELRASTGL GLEHPDTELP LSLAAEPGPA RSPRGSTDLS DQGASPGVLL DPACQKTDHF LSCYPWFHGP ISRVRAAQLV QLQGPDAHGV FLVRQSESRR GEYVLTFLNLQ GRAKHLRLVL TERGQCRVQH LHFPSVVDML RHFQRSPIPL ECGAACDVRL SGYVVVLSQA PGSSNTVLFPSLPHWDSSEL GPHLSSVGC PPSHGAEALP GQVTPPEQIF HLVPSPEELA NSLRQLELES VSSARDSDYD MDSSSRGHLR AIDNQYTPLS QLCREADV
Source	Mammalian cell
Target Names	Sh2b3
Protein Names	Recommended name: SH2B adapter protein 3 Alternative name(s): Lymphocyte adapter protein Lymphocyte-specific adapter protein Lnk Signal transduction protein Lnk
Expression Region	1-548
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.