



# Recombinant DNA polymerase III subunit beta (dnaN)

<b>Product Code</b>	CSB-EP338470EYY-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P22838
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Proteus mirabilis
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKFIIEREQL LKPLQQVSGP LGGRPTLPIL GNLLKVTEN TLSLTGTDLE MEMMARVLSL QSHEIGATTV PARKFFDIWR GLPEGAEISV ELDGDRLLVR SGRSRFSLSL LPASDFPNLD DWQSEVEFTL PQATLKRLIE STQFSMAHQD VRYYLNGMLF ETENTELRTV ATDGHRLAVC AMDIGQSLPG HSVIVPRKGV IELMRLLDGS GESLLQLQIG SNNLRAHVGD FIFTSKLV DG RFPDYRRVLP KNPTKTVIAG CDILKQAFSR AAILSNEKFR GVRINLTNGQ LKITANNPEQ EEAEEIVDVQ YQGEEMEIGF NVSYLLDVLN TLKCEEVKLL LTDAVSSVQV ENVASAAAAY VVMPMRL
<b>Source</b>	E.coli
<b>Target Names</b>	dnaN
<b>Protein Names</b>	Recommended name: DNA polymerase III subunit beta EC= 2.7.7.7
<b>Expression Region</b>	1-367
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.