



Recombinant Transforming protein Maf (V-MAF)

Product Code	CSB-BP338490ASO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23091
Product Type	Recombinant Protein
Immunogen Species	Avian musculoaponeurotic fibrosarcoma virus AS42
Purity	≥85% (SDS-PAGE)
Sequence	MASELAMSGS DLPTSPLAME YVNDFDLMKF EVKKEPVETD RIISQCGRLI AGGSLSTPM STPCSSVPPS PSFSAPSPGS GTDQKTHLED YYWMTGYPQQ LNPEALGFSP EDAVEALINS SHHPLPGAFD GYARGQQLAA AAGGSVPAEE MGSAAAVVSA VIAAAAAQGG APHYHHHHHH PHHGGGGGGG GHPHGAAPGS APPSSASSA AGSGGGGGGG GGGAGGLHHP HHGGGGGGGG LHFDDRFSDE QLVTMSMREL NRQLRGVSKE EVIRLKQKRR TLKNRGYAQS CRFKRVQQRH VLESEKNQLL QQVEHLKQEI SRLVRERDAY KEKYEKLVSN GFRENGSSSD NPSSPEFFMY PRESSTTVM
Source	Baculovirus
Target Names	V-MAF
Protein Names	Recommended name: Transforming protein Maf
Expression Region	1-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.