



Recombinant *Saccharomyces cerevisiae* Glucan 1,3-beta-glucosidase I/II (EXG1)

Product Code	CSB-MP338547SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23776
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	YYDYDHGSLG EPIRGVNIGG WLLLEPYITP SLFEAFRTND DNDEGIPVDE YHFCQYLGKD LAKSRLQSHW STFYQEQDFA NIASQGFNLV RIPIGYWAFQ TLDDDPYVSG LQESYLDQAI GWARNNSLKV WVDLHGAAGS QNGFDNSGLR DSYKFLEDSN LAVTTNVLNY ILKKYSAEEY LDTVIGIELI NEPLGPVLDM DKMKNDYLAP AYEYLRNNIK SDQVIIIHDA FQPYNWDDF MTENDGYWGV TIDHHYQVF ASDQLERSID EHIKVACEWG TGVLNESHWT VCGEFAAALT DCTKWLNSVG FGARYDGSWV NGDQTSSYIG SCANNDIAY WSDERKENTR RYVEAQLDAF EMRGGWIIWC YKTESSLEWD AQRLMFNGLF PQPLTDRKYP NQCGTISN
Source	Mammalian cell
Target Names	EXG1
Protein Names	Recommended name: Glucan 1,3-beta-glucosidase I/II EC= 3.2.1.58 Alternative name(s): Exo-1,3-beta-glucanase I/II Soluble cell wall protein 6
Expression Region	41-448
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.