



# Recombinant Escherichia coli RNA polymerase sigma-54 factor (rpoN)

<b>Product Code</b>	CSB-BP338580ENV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P24255
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKQGLQLRLS QQLAMTPQLQ QAIRLLQLST LELQQELQQA LESNPLLEQI DTHEEIDTRE TQDSETLDTA DALEQKEMPE ELPLDASWDT IYTAGTPSGT SGDYIDDELP VYQGETTQTL QDYLMWQVEL TPFSDTDRAI ATSIVDAVDE TGylTVPLED ILESIGDEEI DIDEVEAVLK RIQRFDVPVGV AAKDLRDCLL IQLSQFDKTT PWLEEARLII SDHLDLLANH DFRTLmRVTR LKEDVLKEAV NLIQSLDRP GQSIQTGEPE YVIPDVLVRK HNGHWTVELN SDSIPRLQIN QHYASMCNNA RNDGDSQFIR SNLQDAKWLI KSLESRNDTL LRVSRCIVEQ QQAFFEQGE EYMKPMVLADI AQAVEMHEST ISRVTTQKYL HSPRGIFELK YFFSSHVNT EGGGEASSTAI RALVKKLIAA ENPAKPLSDS KLTSLLSEQG IMVARRTVAK YRESLSIPPS NQRKQLV
<b>Source</b>	Baculovirus
<b>Target Names</b>	rpoN
<b>Protein Names</b>	Recommended name: RNA polymerase sigma-54 factor
<b>Expression Region</b>	1-477
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.