



Recombinant *Saccharomyces cerevisiae* RNA-binding protein SRO9 (SRO9)

Product Code	CSB-BP338678SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P25567
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSAETAAANT ATAPVPEVQE QESSKSKQVN LTPAPLPTSS PWKLAPTEIP VSTISIEDLD ATRKKKNRTP TPKSSTATKW VPIKASITVS GTKRSGSKNG ASNGNSNKSK NNKTAASSTS SSNANRKKKH HQHNAKKQQQ MKKDGFEAV GEEDSKDATS QENGQSTQQQ QPPHHRNHHH SHHHNSNGPQ RRFHNSNNA GMPQNQGFPF QFKPYQGRNA RNNNNRSKY HNHFHNNQQH PQQPMVKLQQ QFYVPVQVLM AINNIARQIE YYFSEENLTV DNYLRSLKSK DGFAPLSLIS KFYRVVNMSF GGDTNLILAA LREIVANAA TVNVAEGTLA AKEGDNVTGE AKEPSPLDKY FVRSKSWSNW LPETFETEIN IEKELVGDAL DQFMISLPPV PQQEEESSTE LASQEQETKE DSAPVAAGES ESSL
Source	Baculovirus
Target Names	SRO9
Protein Names	Recommended name: RNA-binding protein SRO9 Alternative name(s): Suppressor of RHO3 protein 9
Expression Region	1-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.