



Recombinant Simian foamy virus type 3 Protein Bel-2 (bel2)

Product Code	CSB-BP338835SGQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P27403
Product Type	Recombinant Protein
Immunogen Species	Simian foamy virus type 3 (strain LK3) (SFVagm) (SFV-3)
Purity	>85% (SDS-PAGE)
Sequence	MRDVCLPFLA RENLSNPESG LILLEDTERS HSSLRIGQNA PDGVWPLGNS PILPVVTPWP LCQDHAAPSI WTLDDAYWRG YQDQNLEPPK WLWLCLEDPS GNKYTGQTFL VPPLGLVKIR LYQNLTVVYI CQSIDPWENE NPTGGRRDPT RRYGCRIACD PVYCVKIVWE GNLWDKDDQP CWLIRLKEGH NHGAKELSQR DIKILGESRP YPYGLIGQCP KLQYAIQVKM RVDKAPLTAK VLAVKALHFH RWNICQRENP GIGEGYFSPG YTQALKAYGP QHGSAEQRVW LISTKIVGPQ EKDYWRDAYR WGYFPLVPNK HHPGWTRHLT KFKISRFATP ADIQKIVDEL LPRGASIVTA DGNRYESTRK VHLVNEGTLE EYQAKIRK
Source	Baculovirus
Target Names	bel2
Protein Names	Recommended name: Protein Bel-2
Expression Region	1-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.