



Recombinant Rice tungro bacilliform virus Protein P2

Product Code	CSB-MP338844RER
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P27499
Product Type	Recombinant Protein
Immunogen Species	Rice tungro bacilliform virus (isolate Philippines) (RTBV)
Purity	>85% (SDS-PAGE)
Sequence	MSADYPTFKE ALEKFKNLES DTAGKDKFNW VFTLENIKSA ADVNLASKGL VQLYALQEID KKINNLTTQV SKLPTTSGSS SAGAIVPAGS NTQGQYKAPP KKGIKRKYP A
Source	Mammalian cell
Protein Names	Recommended name: Protein P2 Alternative name(s): ORF 2 P12
Expression Region	1-110
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.