



Recombinant Enterococcus faecalis Mannitol-specific phosphotransferase enzyme IIA component (mtIF)

Product Code	CSB-BP338848ELW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P27547
Product Type	Recombinant Protein
Immunogen Species	Enterococcus faecalis (strain ATCC 700802 / V583)
Purity	>85% (SDS-PAGE)
Sequence	MENLTNISIE LNQQFNTKEE AIRFCGEKLV EAGCVEPAYI EAMIERDQLL SVYMGNFIAI PHGTEEAKKL VKKSGICVVQ VPEGVNFGTE EDEKIATVLF GIAGVGEEHL QLVQQIALYC SDMDNVVQLA DALSKEEITE NLAIA
Source	Baculovirus
Target Names	mtIF
Protein Names	Recommended name: Mannitol-specific phosphotransferase enzyme IIA component EC= 2.7.1.- Alternative name(s): EIIA-Mtl EIII-Mtl PTS system mannitol-specific EIIA component
Expression Region	1-145
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.