



Recombinant Uncharacterized 31.7 kDa protein in traX-finO intergenic region

Product Code	CSB-MP339013ENL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P29368
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	≥85% (SDS-PAGE)
Sequence	MKITDHLKSE GIALTFRVPE GNIKHPLIIL CHGFCGIRNV LLPCFANAFT EAGFATITFD YRGFGESDGE RGRLVPAMQT EDIISVINWA EKQECIDNQR IGLWGTSLGG GHVFSARAQD QRVKCVSQL AFADGDVLVT GEMNESERAS FLSTLNKMAE KKKNTGKEMF VGVTRVLSDN ESKVFFEKVK GQYPEMDIKI PFLTVMETLQ YKPAESAANKV QCPVLIVIAQ QDSVNPPEQG KALYDAVASG TKELYEEADA CHYDIYEGAF FERVAAVQTQ WFKKHL
Source	Mammalian cell
Protein Names	Recommended name: Uncharacterized 31.7 kDa protein in traX-finO intergenic region
Expression Region	1-286
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.