



Recombinant Escherichia coli 2-iminoacetate synthase (thiH)

Product Code	CSB-YP339070ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P30140
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MKTFSDRWRQ LDWDDIRLRI NGKTAADVER ALNASQLTRD DMMALLSPAA SGYLEQLAQR AQLRTRQRFQ NTVSFYVPLY LSNLCANDCT YCGFSMSNRI KRKTLDEADI ARESAAIREM GFEHLLLVGT EHQAKVGM DY FRRHLPALRE QFSSLQMEVQ PLAETEY AEL KQLGLDGVMV YQETYHEATY ARHHLKGGKQ DFFWRLETPD RLGRAGIDKI GLGALIGLSD NWRVDSYMVA EHELLWLQQHY WQSRYSVSFP RLRPCTGGIE PASIMDERQL VQTICAFRL APEIELSLST RESPWFRDRV IPLAINN VSA FSKTQPGGYA DNHPELEQFS PHDDRRPEAV AAALTAQGLQ PVWKDWDSYL GRASQRL
Source	Yeast
Target Names	thiH
Protein Names	Recommended name: 2-iminoacetate synthase EC= 4.1.99.19 Alternative name(s): Dehydroglycine synthase Tyrosine lyase
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.