



Recombinant Drosophila melanogaster Fork head domain-containing protein FD4 (fd96Ca)

Product Code	CSB-BP339227DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32028
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MPRPSRESYG EQKPPYSYIS LTAMAIWSSP EKMLPLSDIY KFITDRFPYY RKNTQRWQNS LRHNLDFNDC FIKVPRRPDR PGKGAYWALH PQAFDMFENG SLLRRRKRFK LHKNDKDLLN EELTALANLN RFFFTTRNGG SAAHMSPLDM NNAAMRLDP LPRSTAHMPN SLGPGVPLPH VMPASMSGAD HTNLADMGLT NLPALTSSEI EGPLSLRPKR SFTIESLITP DKPEHPSEDE DDEDDRVDID VVECSGISRY PTPAASEEY MSASRSSRTE DPLPPMHTIN AGAHVPFLHY ATGANVAGLP ASGIPNSPTT YELAISHPLF MMAAPIANMH NIYYNNVTLV APAQQYRSPE VQNRIDNDR TI
Source	Baculovirus
Target Names	fd96Ca
Protein Names	Recommended name: Fork head domain-containing protein FD4
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.