



Recombinant *Saccharomyces cerevisiae* Pre-mRNA-splicing factor PRP21 (PRP21)

Product Code	CSB-MP339270SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32524
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MEP EDTQLKE DIKTTVNYIK QHGVEFENKL LEDERFSFIK KDDPLHEYTYT KLMNEPTDTV SGEDNDRKSE REIARPPDFL FSQYDTGISR RDMEVIKLT RYYAKDKSIV EQMISKDGEA RLNFMNSSHP LHKTFTDFVA QYKRVSFTG QEIKKSKRTI LDNCFERTQY WEFEKDKDRE HDKLVELCKI QFAAIPWDKF TQVAKFSIPE DTEIFEGSLD LEQMRLRRVQ TGIKLFDSIK PTNEEEKIVS DQGKQKGGDS KGKKRKIRAV GETRLKSKK
Source	Mammalian cell
Target Names	PRP21
Protein Names	Recommended name: Pre-mRNA-splicing factor PRP21
Expression Region	1-280
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.