



Recombinant *Saccharomyces cerevisiae* Nuclear and cytoplasmic polyadenylated RNA-binding protein PUB1 (PUB1)

Product Code	CSB-MP339279SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32588
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	SENNEEQHQ QQQQQQPVAV ETPSAVEAPA SADPSSEQSV AVEGNSEQAE DNQGENDPSV VPANAITGGR ETSRVLVYG NLDKAITEDI LKQYFQVGGP IANIKIMIDK NNKNVNYAFV EYHQSHDANI ALQTLNGKQI ENNIVKINWA FQSQQSSDD TFNLFVGD LN VNVDDETLRN AFKDFPSYLS GHVMWDMQTG SSRGYGFVSF TSQDDAQNAM DSMQGQDLNG RPLRINWAAK RDNNNNNNYQ QRRNYGNNNR GGFRQYNSNN NNNMNMGMNM NMNMNMNNSR GMPPSSMGMP IGAMPLPSQG QPQQSQTIGL PPQVNPQAVD HIIRSAPPRV TTAYIGNIPH FATEADLIPL FQNFILDF KHYPEKGCCF IKYDTHEQAA VCIVALANFP FQGRNLRRTGW GKERSNFMPQ QQQQGGQPLI MNDQQQPVMS EQQQQQQQQQ QQQ
Source	Mammalian cell
Target Names	PUB1
Protein Names	Recommended name: Nuclear and cytoplasmic polyadenylated RNA-binding protein PUB1 Alternative name(s): ARS consensus-binding protein ACBP-60 Poly uridylyate-binding protein Short name= Poly(U)-binding protein
Expression Region	2-453
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.