



Recombinant *Saccharomyces cerevisiae* 54S ribosomal protein RML2, mitochondrial (RML2)

Product Code	CSB-YP339282SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32611
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	QENTSKD DSSLNITPRL LKIIPNDTDI VTLEKQDELI KRRRKLSKEV TQMKRLKPVS PGLRWYRSPI YPYLYKGRPV RALTVVRKKH GGRNNSGKIT VRHQGGGHRN RTRLIDFNW EGGAQTVQRI EYDPGRSSHI ALLKHNTTGE LSYIIACDGL RPGDVVESFR RGIPQTLLE MGGKVDPAIL SVKTTQRGNC LPISMIPIGT IIHNVGITPV GPGKFCRSAG TYARVLAKLP EKKKAIVRLQ SGEHRYVSLE AVATIGVVSNI DHQNRSLGK AGRSRWLGIPT VTRGVAMNK CDHPHGGGRG KSKSNKLSMS PWGQLAKGYK TRRGKNQNRM KVKDRPRGKD ARL
Source	Yeast
Target Names	RML2
Protein Names	Recommended name: 54S ribosomal protein RML2, mitochondrial Short name=L2
Expression Region	44-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.