



Recombinant Human Zinc finger protein 193 (ZNF193)

Product Code	CSB-BP026569HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O15535
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNTNSKEVLS LGVQVPEAWE ELLTMKVEAK SHLQWQESRL KRSNPLAREI FRRHFRQLCY QETPGPREAL TRLQELCYQW LRPHVSTKEQ ILDLLVLEQF LSILPKELQG WVREHCPESE EEAILLEDL ERELDEPQHE MVAHRHRQEV LCKEMVPLAE QTPLTLQSQP KEPQLTCDSDA QKCHSIGETD EVTKTEDREL VLRKDCPKIV EPHGKMFNEQ TWEVSQQDPS HGEVGEHKDR IERQWGNLLG EGQHKCDECG KSFTQSSGLI RHQRIHTGER PYECNECGKA FSRSSGLFNH RGIHNIQKRY HCKEKGKVFV QSAGLIQHQR IHKGEKPYQC SQCSKSYSRR SFLIEHQRSH TGERPHQCIE CGKSFNRHCN LIRHQKIHTV AELV
Source	Baculovirus
Target Names	ZSCAN9
Protein Names	Recommended name: Zinc finger protein 193 Alternative name(s): Cell proliferation-inducing gene 12 protein PRD51 Zinc finger and SCAN domain-containing protein 9
Expression Region	1-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.