



Recombinant *Saccharomyces cerevisiae* GTP-binding nuclear protein GSP2/CNR2 (GSP2)

Product Code	CSB-EP339300SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32836
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	SAPAQNNAE VPTFKLVLVG DGGTGKTTFV KRHLTGEFEK KYIATIGVEV HPLSFYTNFG EIKFDVWDTA GQEKFGGLRD GYYINAQCAI IMFVTSRIT YKNVFNWHRD LVRVCENIPI VLCGNKVDVK ERKVKAKTIT FHRKKNLQYY DISAKSNYNF EKPFLWLARK LAGNPQLEFV ASPALAPPEV QVDEQLMHQY QQEMDQATAL PLPDEDDADL
Source	E.coli
Target Names	GSP2
Protein Names	Recommended name: GTP-binding nuclear protein GSP2/CNR2
Expression Region	2-220
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.