



Recombinant *Saccharomyces cerevisiae* 26S protease regulatory subunit 6B homolog (RPT3)

Product Code	CSB-EP339339SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P33298
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MEELGIVTPV EKAVEEKPAV KSYASLLAQL NGTVNNNSAL SNVNSDIYFK LKKLEKEYEL LTLQEDYIKD EQRHLKRELK RAQEEVKRIQ SVPLVIGQFL EPIDQNTGIV SSTTGMSYVV RILSTLDREL LKPSMSVALH RHSNALVDIL PPDSOSSISV MGENEKPDVT YADVGGLDMQ KQEIREAVEL PLVQADLYEQ IGIDPPRGVL LYGPPGTGKT MLVKAVANST KAAFIRVNGS EFVHKYLGEG PRMVRDVFRL ARENAPSIIF IDEVDSIATK RFDAQTGS DR EVQRILIELL TQMDGFDQST NVKVIMATNR ADTLDPALLR PGRLDRKIEF PSLRDRRERR LIFGTIASKM SLAPEADLDS LIIRNDSLGS AVIAAIMQEA GLRAVRKNRY VILQSDLEEA YATQVKTDNT VDKFDFYK
Source	E.coli
Target Names	RPT3
Protein Names	Recommended name: 26S protease regulatory subunit 6B homolog Alternative name(s): Protein YNT1 Tat-binding homolog 2
Expression Region	1-428
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.