



Recombinant Escherichia coli Uncharacterized protein yehM (yehM)

Product Code	CSB-MP339344ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P33349
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSEPLIVGIR HHSPACARLV KSLIESQRPR YVLIEGPADF NDRVDELFLA HQLPVAIYSY CQYQDGAAPG RGAWTPFAEF SPEWQALQAA RRIQAQTYFI DLPCWAQSEE EDDSPDTQDE SQALLLRATR MDNSDTLWDH LFEDESQQTA LPSALAHYFA QLRGDASGDA LNRQREAFMA RWIGWAMQQN NGDVLVCCGG WHAPALAKMW RECPQKINKP ELPSLADAVT GCYLTPYSEK RLDVLGAGYLS GMPAPVWQNW CWQWGLQKAG EQLLKTILTR LRQHKLPAST ADMAAAHLHA MALAQLRGHT LPLRTDWLDA IAGSLIKEAL NAPLPWSYRG VIHPDTPIL LTLIDTLAGD GFGKLAPSTP QPPLPKDVT C ELERTAISLP AELTLNRFTP DGLAQSQVLH RLAILEIPGI VRQQGSTLTL AGNGEERWKL TRPLSQHAAL IEAACFGATL QEAARNKLEA DMLDAGGIGS ITTCLSQAAL AGLASFSQQL LEQLTLIAQ ENQFAEMGQA LEVLYALWRL DEISGMQGAQ ILQTTLCATI DRTLWLCE SN GRPDEKEFHA HLHSWQALCH ILRDLHSGVN LPGVSLSAAV ALLERRSQAI HAPALDRGAA LGALMRLEHP NASAEAALTM LAQLSPAQSG EALHGLLALA RHQLACQPAF IAGFSSHLNQ LSEADFINAL PDLRAAMAWL PPRERGT LAH QVLEHYQLAQ LPVSALQMPL HCPPQAI AHH QQLEQQALAS LQNWGVFHV
Source	Mammalian cell
Target Names	yehM
Protein Names	Recommended name: Uncharacterized protein yehM
Expression Region	1-759
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.