



# Recombinant *Saccharomyces cerevisiae* Protein phosphatase 2C homolog 3 (PTC3)

<b>Product Code</b>	CSB-EP339428SVG-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P34221
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGQILSNPII DKEHHS GTDC LTAFG L CAMQ GWRMSMEDAH IVEPNLLAES DEEH LAFYGI FDGHGGSSVA EFCGSKMISI LKKQESFKSG MLEQCLIDTF LATDV ELLKD EKLKDDHSGC TATVILVSQL KKLICANSG DSRTVLSTGG NSKAMSF D HK PTL LSEK SRI VAADGFVEMD RVNGNLALSR AIGDFEFKSN TKLGPHEQVV TCV PDIICHN LNYDEDEFVI LACDGIWDCL TSQECVDLVH YGISQGNMTL SDISSRIVDV CCSPTTEGSG IGC DNMSISI VALLKENESE SQWFERMR SK NYNIQT SFVQ RRKSIFDFHD FSDDDNEVFA ITTKKLQDRL NRSKDND DME IDDLDELGS SATPSKLSGE DRTGPIDLFS LEALLEAGIQ IRQRPS DSD GNTSYFH GAS LSDMLASLSN AAAGETEPND ADDNDDNDGE ENGNENAKK GSKIEEIE
<b>Source</b>	E.coli
<b>Target Names</b>	PTC3
<b>Protein Names</b>	Recommended name: Protein phosphatase 2C homolog 3 Short name= PP2C-3 EC= 3.1.3.16
<b>Expression Region</b>	1-468
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.