



# Recombinant *Saccharomyces cerevisiae* DNA replication regulator SLD2 (SLD2)

<b>Product Code</b>	CSB-MP339432SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P34252
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MYSFELDKLK IELKTWEHDF IDKNKREPTR DDIKSLRTVR QMYKQYSTLK KKQSLQRQKV DTQESVELPA HKKDHDEVVE IGPTPQVYVK AISIFDMNLS PIKPIYMTFT NNIDVNNDNS KTISNESSPR KTILLKSSPA DRTLVAEPIS SVKRQLNFQM LNASSTRPT SSPCKNRNGK LVEIKKCSPT INPPLESGKP SGYYGPN SPL KLDEENIHLN ISLNSSTKRR LQIAYPSLQK TPSKDQADIS TSFSPSLIR RPLTKSLIEL AREHTEIVKE FGV LQEEDIE EEEEEEGEEN GYDEKNHEDD FGLEDELIRP KVVKDIFQED DDNDDSQARE DTFIRKRPKR RKVIRRLRDN DPETETAGFE RDVHKELVKL KRRKVAEFLG STSQISDTEF EHDPEASSGV VSSEQKPTAK RKGRKKYNLV SNNFRRLKLP KKNRFSNGRW GRR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SLD2
<b>Protein Names</b>	Recommended name: DNA replication regulator SLD2 Alternative name(s): DNA replication and checkpoint protein 1
<b>Expression Region</b>	1-453
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.