



Recombinant Uncharacterized protein C30A5.4 (C30A5.4)

Product Code	CSB-EP339445CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34350
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSPMNRQRKN KSNVLNEKDE RPGCIIKEKS SVRSMAKTQP TSGRTISLSE INKPIALSJK NTEIKAAVEI QNRLTKKMAK NSERRKKMAH GPEMTAQNRT QERIDFVEFL DRHYQVMQAG VSTVKVHKMH YTKEEFENVI YEAQTIFSSE KALVDIDPPC VVVGNLHGQF NDLINMFILL GRPPETVYVF TECVTSIPKC GEEIWALFQR CFNNLPISAL IATKILCMHG GLSPALTCLD ELRNHPKPIR NPFRGIVNDM LWADPDPSVF EWKTSSRGSF FTFGTNVIDD VCKRLGVELI IRAHQMCFDG YWVLSGRKLI TIFSAPMYCN FYKNTGCVLK VDETLGIQTV AFVPESENIE KIIIEEMNRVW DISIDCIE
Source	E.coli
Target Names	C30A5.4
Protein Names	Recommended name: Uncharacterized protein C30A5.4
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.