



Recombinant Putative zinc finger protein F44E2.7 (F44E2.7)

Product Code	CSB-BP339456CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34437
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MVKPTSSLLQ RYKNNVRVDP STLNHLYLCYY CGKTLSDRLE YQQHMLKVHE VMSQSYQLWS DHTMDTLGLI WSESKRSYDR DTSCSLSVPA SPMSRESRNR NSDDDDYYDYV YEDEKPKRRR IDVEQSRRSV AAAAEKHRKL LEEQRRKAEA EYERKRKADE KKKEKAIKDE EKRKAEDLRK CLQLQQQKTA AAAQIRETAE KAAVARRMST FEVGESSEQL AKPEGEKRKR RRTESRWREI ESDSDKPEVT ALVKKILEES KKKEASSQEV EDADLVPELS TRKPYDNTEH VSKLCKVCKK GPHYTFANLF EHYQDLHNAT IKSLHYYGFN GNKLIGKKNL IERDHCQRCV IKFPRARDYF AHMIKHHVYE SVRCQLDFEN ATNADVARM MFRDRILTLG YNFKFEQVAD PNLVSDVLEP GQEPSTSAEQ EDPSSLKIVK LEEPELEEQK ALKPESEAME QQEDVHSLNP E
Source	Baculovirus
Target Names	F44E2.7
Protein Names	Recommended name: Putative zinc finger protein F44E2.7
Expression Region	1-471
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.